


INSTRUCTIONS #11 Finger Knuckle Bender Splint

Applying splint with only one end of the rubber band attached to hook (no dynamic force to flex P.I.P. joint) will facilitate proper positioning of splint on finger. Slide finger into splint, position wider finger band on top of finger and behind D.I.P. joint. Position bar under P.I.P. joint. Complete attaching rubber bands. Evaluate fit, function, and instruct patient on use. The splint is designed to dynamically flex finger while allowing active extension. Rubber bands can be added or subtracted and the wires and bands can be bent to modify fit and function. Remove splint from patient to modify.

Disconnecting rubber bands will make splint removal easier and prepare splint for reapplication.

Clean with mild soap and damp towel. Air dry completely before reapplying. Immediately contact your health care professional if you are experiencing loss of sensation, swelling, skin irritation, or have any questions or concerns.



Bunnell #11 Finger Knuckle Bender to flex proximal interphalangeal joint

CAUTION: Federal Law restricts this device to sale by or on the order of a physician



Single patient use